



HRV Coach Pro v2.1

Autonomous Clinical Intelligence

Khan Saquib , Liu TzuEn Andrew , Peng ToChen
National Cheng Kung University



Project Overview

- **HRV Coach Pro v2.1:** An autonomous agentic system for clinical HRV analysis.
- **Core Function:** Converts raw noisy ECG → validated HRV metrics → clinical-style insights.
- **Agentic Approach:** Utilizes a self-correcting workflow to adapt to signal quality.
- **Goal:** Provide scalable, reliable HRV analysis without requiring expert manual intervention.





The Signal-to-Insight Gap

- **The Pipeline:** ECG \rightarrow R-peaks detection \rightarrow RR Intervals \rightarrow HRV Metrics.
- **The Challenges:** Noise, baseline drift, motion artifacts, and EMG interference often cause detection failure.
- **The Risk:** Incorrect peak detection leads to invalid RR intervals and distorted HRV results.
- **Interpretation:** Standard metrics (RMSSD, SDNN) are difficult for non-experts to understand without context.

Agentic Workflow & Architecture



NORDLING LAB



Sense → Decide → Act → Verify

- **Data Ingestion:** PhysioNet + CSV support.
- **Signal Processing:** Filtering, peak detection, and RR calculation.
- **Strategy Control:** Validation-driven feedback loop.

Hybrid Intelligence

- **Deterministic:** NeuroKit2 & SciPy for signal processing.
- **Probabilistic:** DeepSeek V3.2 for metric interpretation.
- **Safety:** Separation of computation vs. interpretation.



Strategy Pattern: Why Multiple Pipelines?

- **Adaptability:** A single pipeline fails across different signal conditions.
- **Strategy Library:**
 - **Strategy A:** Standard filtering for clean clinical-grade ECG.
 - **Strategy B:** Aggressive filtering for noisy wearable data.
 - **Strategy D:** Localized handling for 50 Hz power-line interference.
- **Selection:** The agent tests strategies and selects the best one based on a validation grade.



Decision Loop & Validation Logic

- 1 Sensing:** Checks for baseline wander, high-frequency noise, and sampling consistency.
- 2 Decision:** Try Strategy A → Detect Peaks → Validate results.
- 3 Backtracking:** If Grade is "C" or "Reject", switch to Strategy B or D.
- 4 Criteria:** Physiological HR range, RR consistency, and outlier artifact rejection.
- 5 Output:** Only accept results that pass physiological plausibility.



Challenges & Results

- **NumPy 2.0 Compatibility:** Fixed breaking changes in NeuroKit2 by implementing a compatibility shim for trapezoid.
- **Local Data:** Resolved 50 Hz noise issues via Strategy D and timestamp-based sampling inference.
- **Key Deliverables:**
 - Clean ECG plots with labeled R-peaks.
 - Signal quality grading and HRV summary.
 - Downloadable PDF reports with LLM-generated interpretations.



Future Work

- Integration of local LLM models to enhance data privacy.
- Real-time edge deployment for continuous health monitoring.
- Expansion of the strategy library to support more diverse sensor types.



Team Contribution

Member	Contribution	Responsibilities
Khan	40%	Code Development, Video Recording
Liu	30%	Documentation, Hosting Meeting
Peng	30%	Video Editing, Slide Creation

